

# Endo Optiks® Ophthalmic Endoscopes

See. Diagnose. Treat.



Endo Optiks® ophthalmic endoscopes and probes offer a bright new perspective for treatment needs across many surgical specialties. Beaver-Visitec is proud to offer ophthalmic microendoscopy as a powerful surgical solution for minimally invasive glaucoma surgery (MIGS), combined phaco-ECP and vitreoretinal surgery.

Endoscopes and probes can be used with the E2 Ophthalmic Laser and Endoscopy System and the E4 Ophthalmic Endoscopy System. Endoscopes and Probes are available in a variety of sizes and resolutions to meet all of your surgical needs.

Key features include:

- All endoscopes and probes are autoclavable.
- Sizes range from 25 gauge to 18 gauge.
- High resolution image display.
- Ultra wide field imaging.
- Ultra wide field illumination.
- Available in straight or curved tip configurations.

# Ophthalmic Endoscopes and Probes

Our patented ophthalmic laser endoscopes and probes are available in triple function, dual function and single function configurations, delivering a combination of visualized laser application, high resolution imaging and wide field illumination, through one single microprobe.

## Triple Function Ophthalmic Laser Endoscopes

Guage	Cannula	Resolution	Laser	Camera	Light	Custom Laser Connector	E2 Compatibility	E4 Compatibility	Item #
23	Straight	High (10k pixel)	✓	✓	✓	✓	✓	✓	OME230SMAHRXA
23	Straight	Standard (6k pixel)	✓	✓	✓	✓	✓	✓	OME230SMAXA
20	Straight	High (10K pixel)	✓	✓	✓	—	✓	—	OME200SMAEA
20	Straight	High (10k pixel)	✓	✓	✓	✓	—	✓	OME200SMAXA
19	Straight	High (10k pixel)	✓	✓	✓	—	✓	—	OME200SMAHRA
19	Curved	Ultra High (17k pixel)	✓	✓	✓	—	✓	—	OME200SMAHRVA
19	Straight	Ultra High (17k pixel)	✓	✓	✓	✓	—	✓	OME200SMAHRXA
19	Curved	Ultra High (17k pixel)	✓	✓	✓	✓	—	✓	OME200SMAHRXVA
19	Curved	High (10k pixel)	✓	✓	✓	—	✓	—	OME200SMAVA
19	Curved	High (10k pixel)	✓	✓	✓	✓	—	✓	OME200SMAXVA

## Dual Function Ophthalmic Endoscopes

Guage	Cannula	Resolution	Laser	Camera	Light	Custom Laser Connector	E2 Compatibility	E4 Compatibility	Item #
23	Straight	Standard (6k pixel)	—	✓	✓	—	✓	✓	OME230LA
23	Straight	High (10k pixel)	—	✓	✓	—	✓	✓	OME230HRA
20	Straight	High (10K pixel)	—	✓	✓	—	✓	✓	OME200LA
20	Straight	Ultra High (17k pixel)	—	✓	✓	—	✓	✓	OME200HRA
19	Curved	High (10k pixel)	—	✓	✓	—	✓	✓	OME200LVA

## Dual Function Ophthalmic Laser Probes

Guage	Cannula	Resolution	Laser	Camera	Light	Custom Laser Connector	E2 Compatibility	E4 Compatibility	Item #
23	Straight	—	✓	—	✓	—	✓	✓	VPH230LA
20	Straight	—	✓	—	✓	—	✓	✓	OME20LA

## Single Function Ophthalmic Probes

Guage	Cannula	Resolution	Laser	Camera	Light	Custom Laser Connector	E2 Compatibility	E4 Compatibility	Item #
25	Straight	—	✓	—	—	—	✓	—	VPH250A
23	Straight	—	✓	—	—	—	✓	—	VPH230A
23	Straight	—	—	—	Wide Field Illumination	—	✓	✓	VP230WA
20	Straight	—	—	—	Wide Field Illumination	—	✓	✓	VP200WA



### Beaver-Visitec International, Inc.

411 Waverley Oaks Road, Waltham, MA 02452  
 Tel: 1.866.906.8080 Fax: 1.866.906.4304  
 www.beaver-visitec.com

### Beaver-Visitec International, Sales Limited

85c Park Drive, Milton Park,  
 Abingdon, Oxfordshire, OX14 4RY, UK  
 Tel: 44.1865.601.256 Fax: 44.1865.595.761