


Parasol® Punctal Occluder

Visibly Different

One glance at the Parasol® Punctal Occluder and you'll recognize a breakthrough product offering distinct benefits for you and your dry eye patients:



Easy To Insert

The Parasol's hollow nose collapses on insertion, eliminating the need for punctal dilation in most cases. It's easier on you and more comfortable for your patients.

Simple Sizing

Two sizes, small and medium, fit virtually all patients. X-small and large sizes are available for uncommonly sized puncta.

Designed for Excellent Retention*

Size Comparison Chart

| Traditional Plugs | 0.2 mm | 0.3 mm | 0.4 mm | 0.5 mm | 0.6 mm | 0.7 mm | 0.8 mm | 0.9 mm + |
|-------------------|---------|--------|--------|--------|--------|--------|--------|----------|
| Parasol | X-SMALL | SMALL | | | MEDIUM | | | LARGE |

Small and medium sizes should accommodate 95% of all patients. BVI recommends undersizing when in doubt.

| Parasol® Sterile Pre-Loaded Punctal Occluders | Quantity | Item # |
|---|-------------------------------|-----------|
| X-Small (0.25mm-0.35mm) | 1 Pair | A14-200 |
| Small (0.35mm-0.65mm) | 1 Pair | A14-201 |
| Medium (0.6mm-0.85mm) | 1 Pair | A14-202 |
| Large (0.9mm+) | 1 Pair | A14-203 |
| Parasol® Sterile Quickload™ Punctal Occluders | | |
| Small (0.35mm-0.65mm) | 8 Pairs | A14-201QL |
| Medium (0.6mm-0.85mm) | 8 Pairs | A14-202QL |
| Large (0.9mm+) | 8 Pairs | A14-203QL |
| Combo Pack | 4 Small Pairs, 4 Medium Pairs | A14-205QL |
| Parasol® Non-Sterile Bulk Punctal Occluders | | |
| Small (0.35mm-0.65mm) | 10 Pairs | A14-201BN |
| Medium (0.6mm-0.85mm) | 10 Pairs | A14-202BN |
| Large (0.9mm+) | 10 Pairs | A14-203BN |
| Large (0.9mm+) | 5 Pairs | A20-203BN |
| Combo Pack | 5 Small Pairs, 5 Medium Pairs | A14-209BN |

* McCabe, C. (2009) Punctal occlusion reduces dry eye symptoms and improves vision. Review of Ophthalmology, 16(1), 55-58.

** Certain conditions apply. Call for details. Parasol® U.S. Patent No. 5,830,171

How To Order: Call 866-906-8080 or www.beaver-visitec.com

Beaver-Visitec International, Inc.

411 Waverley Oaks Road, Waltham, MA 02452
Tel: 1.866.906.8080 Fax: 1.866.906.4304
www.beaver-visitec.com

Beaver-Visitec International, Limited

85c Park Drive, Milton Park,
Abingdon, Oxfordshire, OX14 4RY, UK
Tel: 44.1865.601.256 Fax: 44.1865.595.761

