

Visitec® FLACS Cannula




Innovative tip design allows multiple functions to be performed:

1. Opens all femto incisions (side-port, main & LRIs)
2. Allows for injection of OVD and/or other solutions
3. Checks the capsulotomy has been released

Designed to improve:

- Surgical efficiency
- Clinical performance
- Patient safety



Visitec® FLACS Cannula

Excellent Performance


- Thin and blunt leading edge ensures incisions are opened cleanly, without distortion or stretching

Increased Surgical Efficiency


- The multi-functional design permits fewer hand offs between surgeon & scrub nurse, thereby reducing operating time

Patient Safety


- Cannula is single-use to eradicate the risk of cross contamination from reusable instrument









1. Opening of femto incisions — phaco, paracentesis & LRIs



2. Injection of OVD and/or other solutions



3. Check capsulotomy has been fully released from the capsule

FLACS Cannula		Item #	Quantity
FLACS Cannula [Schweitzer] 45 degrees, 9 mm from end, 25G.	 	585200	10 per box
Related Products for Femto Cataract Surgery		Item #	Quantity
Nucleus Hydrodissector [Pearce] 35 degrees, 8 mm from end, 27G.	 	585036	10 per box
Capsulorhexis Forceps with Grasping Tips Forceps have blunt tips with wide platforms which are ideally suited to removing the femto capsulotomy.	 	581486	10 per box

How To Order: Call 866-906-8080 or www.beaver-visitec.com



Beaver-Visitec International, Inc.
 411 Waverley Oaks Road, Waltham, MA 02452
 Tel: 1.866.906.8080 Fax: 1.866.906.4304
www.beaver-visitec.com

Beaver-Visitec International, Sales Limited
 85c Park Drive, Milton Park,
 Abingdon, Oxfordshire, OX14 4RY, UK
 Tel: 44.1865.601.256 Fax: 44.1865.595.761