

# Keracel®

High Density, Lint Free

Keracel products are constructed of a patented high-density formulation of lint free Merocel polyvinyl acetal (PVA) sponge. Designed to benefit both the patient and surgeon, each device is absorbent and firm, allowing gentle manipulation of delicate tissues during LASIK and other corneal procedures.

## Keracel® Spears



Triangle-shaped high density PVA spear set on a malleable polypropylene handle.

Description	Item #
Eye Spear, 60 Count, 3 Spears/pouch, 20/box	403101
Eye Spear, 200 Count, 10 Spears/pouch, 20/box	410101

## Keracel® Lasik Kits



Description	Item #
Contains 5 Keracel Eye Spears, 1 Chayet Drainage Ring, 1 Corneal Shield, 1 Instrument Wipe, 10/box	400131
Contains 3 Keracel Eye Spears, 5 Merocel Eye Spears, 20/box	453101

## Visitec® Spears



Description	Item #
Lasik Expanded Eye Spear, 40 count, 2 spears/pouch, 20/box	581709
Lasik PVA Spear (Blue), 100 count, 10 spears/pouch, 10/box	581751
Lasik PVA Spear (Orange), 100 count, 10 spears/pouch, 10/box	581761

For a complete listing of Fluid Control Products, refer to our product catalog.




**Beaver-Visitec International, Inc.**  
 411 Waverley Oaks Road, Waltham, MA 02452  
 Tel: 1.866.906.8080 Fax: 1.866.906.4304  
[www.beaver-visitec.com](http://www.beaver-visitec.com)

# Fluid Control Products

Unparalleled Fluid Control: Merocel<sup>®</sup>, Weck-Cel<sup>®</sup>, Keracel<sup>®</sup>

Beaver-Visitec International's range of fluid control products are designed to meet the needs of clinicians for a variety of ophthalmic procedures.



For ordering information call:  
1.866.906.8080/1.781.906.8080  
or visit [www.beaver-visitec.com](http://www.beaver-visitec.com)

**bvi** Beaver Visitec  
*Keeping Your Vision in Sight*

# Merocel®

Highly Absorbent, Fast Wicking, Lint Free

Merocel® ophthalmic products are designed to provide the surgeon with superior fiber free products for unparalleled fluid control and safe cleaning of fine microsurgical instruments. The exclusive Merocel material is highly absorbent, fast wicking, lint free, polyvinyl acetal (PVA).

## Merocel® Spears/Sponge Points



Triangle shaped spearhead set on a malleable, polypropylene handle. Spear design allows for use in delicate, surgical areas.

Description	Item #
Eye Spears, 180 count, 10 spears/pouch, 18 pouches/box	400101
Eye Spears, 50 count, 5 spears/pouch, 10 pouches/box	405101
Sponge Points, 400 count, 20 points/pouch, 20 pouches/box	400115
Sponge Points, 200 count, 10 points/pouch, 20 pouches/box	410115

## Merocel® Kits

Merocel® Eye Kits are designed to provide superior fluid control solutions for your entire procedure.

Description	Item #
Contains 1 Wipe and Wick, 1 Corneal Light Shield, 20/box	400117
Contains 1 Wipe, 1 Eye Drain, 1 Corneal Light Shield, 20/box	223622
Contains 5 Eye Spears, 1 Wipe, 1 Corneal Light Shield, 80cc Eye Drain with collection bag, 10/box	400113

## Merocel® Wipes



Description	Item #
Instrument Wipe, 3.625" x 3", 20/box	400200
Instrument Wipe, 3" x 3", 20/box	223620
Instrument Wipe, 2" x 2", Adhesive Back, 20/box	400204
Instrument Wipe, 1.5625" x 3.25", Adhesive Back, 20/box	223625
Instrument Wipe and Wick, 20/box	400210
Instrument Wipe and Eye Drain, 3" x 3", 20/box	223621

## Merocel® Eye Drain with Collection Bag


Description	Item #
6" Eye Wick, with 80cc collection bag. 1 set per packet, 20 packets/box	400108
8" Eye Wick, with 400cc collection bag. 1 set per packet, 10 packets/box	400109

# Weck-Cel®

Fast Absorbency And Affinity For Fluids

The standard in eye surgical fluid control, Weck-Cel® products are made of highly absorbent, natural, sterile cellulose material. The Weck-Cel® cellulose material is designed to provide the surgeon with superior absorption and affinity for fluid, all the while maintaining its rigidity during the wicking process.

## Weck-Cel® Spears



Triangle shaped spearhead, set on a malleable, polypropylene handle, is designed for use in delicate, surgical areas. Designed to maintain rigidity during wicking process.

Description	Item #
Eye Spears, 144 count, 6 spears/pouch, 24 pouches/box	0008680
Eye Spears, 180 count, 10 spears/pouch, 18 pouches/box	0008685

## Weck-Cel® Sponge Points & Strips

Description	Item #
Points, 240 count, 20 points/pouch, 12 pouches/box	0008661
2 Points and 6 Strips/pouch, 24 pouches/box	0008660

## Visitec® Eye Sponge



Description	Item #
Eye Sponge, 100 count, 5 sponges/pouch, 20/box	581089

For a complete listing of Fluid Control Products, refer to our product catalog.